

Test of Probe

PH MASDVT1000S6 BL1.3

N/C/H NO_I/E ,

H170062/0001

1 Introduction

The probe was tested for achievable pulse width/RF field strength, resolution/line shape, signal-to-noise, stability during decoupling and CP, and double CP efficiency and stability. The tests were performed on 18th of August, 2020. This report was finished on the 13th of October, 2020.

2 Test Summary

2.1 Probe Performance

RF field strength and power requirements			
Irradiation type	pulse width/ μ s	RF field strength/kHz	Power (TopSpin) / W
¹ H decoupling	1.8	140	45
¹³ C CP	3.0	83	40
¹⁵ N CP	3.5	70	82
Resolution			
nucleus	line width	Sample/parameters	comment
¹³ C	< 2.0 Hz	Adamantane, ns=16, d1=10s	40 kHz spinning, 20 kHz SPINAL64 decoupling
Signal to noise			
nucleus	S/N	Sample/parameters	comment
¹³ C	94	Glycine (¹³ C and ¹⁵ N in natural abundance), 10 kHz spinning, ns=64, aq=14 ms, 140 kHz decoupling	
¹⁵ N	15	Glycine (¹³ C and ¹⁵ N in natural abundance), 10 kHz spinning, ns=64, aq=29ms, 140 kHz decoupling	
CP stability			
nucleus	result	Sample/parameters	comment
¹⁵ N	ok	Glycine (¹³ C and ¹⁵ N in natural abundance), ns=4, aq=29 ms, 140 kHz decoupling, d1=5s 256 spectra with 4scans each at 5s recycle delay and 5.04 s repetition time	test duration: 1h 25 min

Double CP efficiency			
sample type	efficiency	Sample/ parameters	comment
Full rotor	52%	Glycine-1,2- ¹³ C ₂ - ¹⁵ N, full spinner, 10 kHz spinning, ¹⁵ N- ¹³ C contact: 16 ms, 90% to 100% ramped ¹³ C pulse centered at 25 kHz, square ¹⁵ N pulse at 35 kHz, ¹ H- ¹⁵ N contact: 3 ms 50% to 100% ramped ¹ H pulse centered at 80 kHz, square ¹⁵ N pulse at 70 kHz, 140 kHz decoupling	¹³ C at -1 sideband condition in ¹⁵ N- ¹³ C, contact parameters optimized on sample
	44%	Glycine-1,2- ¹³ C ₂ - ¹⁵ N, full spinner, 10 kHz spinning, ¹⁵ N- ¹³ C contact: 20 ms, 90% to 100% ramped ¹³ C pulse centered at 45 kHz, square ¹⁵ N pulse at 35 kHz, ¹ H- ¹⁵ N contact: 3 ms 50% to 100% ramped ¹ H pulse centered at 80 kHz, square ¹⁵ N pulse at 70 kHz, 140 kHz decoupling	¹³ C at +1 sideband condition in ¹⁵ N- ¹³ C, contact parameters optimized on sample
Double CP stability			
nucleus	result	sample/parameters	comment
¹ H, ¹⁵ N, ¹³ C	ok	Glycine-1,2- ¹³ C ₂ - ¹⁵ N, full spinner, 10 kHz spinning, ¹⁵ N- ¹³ C contact: 20 ms, 90% to 100% ramped ¹³ C pulse centered at 55 kHz, square ¹⁵ N pulse at 32 kHz, ¹ H- ¹⁵ N contact: 3 ms, 50% to 100% ramped ¹ H pulse centered at 80 kHz, square ¹⁵ N pulse at 70 kHz, 140 kHz decoupling 512 spectra with 8 scans each at 5 s recycle delay, total number of scans: 4096	¹³ C at +2 sideband condition total experiment time: 5 h 40 min

Notes:

- 1) The ¹⁵N signal to noise test was run at 10 kHz spinning as defined for a 1 GHz system. As the ¹⁵N figure is rather low compared to the results obtained for comparable 800 MHz probes the test were run at 8 kHz spinning. An S/N of 17 was found (using optimized 50% to 100% or 70% to 100% ramped ¹H contact pulses, respectively). This figure is not significantly different from the figure of 15 determined at 10 kHz spinning.
- 2) All double CP efficiencies are determined by comparing the intensity of the C_α signal in the double cross polarization experiment (¹H → ¹⁵N → ¹³C) to the intensity in the single cross polarization experiment (¹H → ¹³C) for the same sample. The single cross polarization experiment was optimized using the standard 50% to 100% ramp shaped ¹H contact.
- 3) The double CP stability test was run at the +2 sideband condition (which is not the most efficient one) for ¹³C deliberately to stress the probe.

3 Experimental

3.1 Spectrometer setup

Instrument: Avance Neo 1000 SB

TopSpin: 4.0.9 (of 2020-05-25 14:09:09)

3.2 Standard Temperature settings

- T(set)= 20°C
- VT=500 l/h
- p (bearing) and p (drive) as set by the MAS III unit (rotation profile: Generic 0.7 mm)
- Regulation on TC2 (probe has a single thermocouple, only)
- Chiller=BCU II, power = strong

3.3 External Filters

Filters were used on all probe channels throughout:

^{15}N channel: ^{15}N LRP filter

^{13}C channel: ^{13}C - ^{23}Na BP filter

^1H channel: ^1H BP filter

3.4 Double Cross Polarization Experiment

This experiment was optimized using a 90% to 100% ramped ^{13}C contact pulse for the $^{15}\text{N} \rightarrow ^{13}\text{C}$ contact, no other shapes (square, tangential) were tried.

The determined maximum efficiencies (-1, +1, and +2 sideband conditions) are good. Long contact pulses were possible indicating that the probe's tuning and matching are stable during the pulses.

4 Figures




Figure 1: Magic angle adjustment: ^{79}Br spectrum of KBr at 25 kHz spinning.




Figure 2: ^{13}C resolution test: ^{13}C spectrum (direct excitation) of adamantane at 40 kHz spinning.

¹³C CP and sensitivity
 glycine (natural abundance), 10 kHz spinning
¹³C CP at 83 kHz (3.0 us) / 2 ms, plw1=40 W
 decoupling: SPINAL64 at 140 kHz (1.8 us) / 14 ms, plw12=45 W
 temperature control on, T(set)= 293 K, VT=500 l/h
 S/N: 94 (ns=64, sinocal, 20 ppm noise)




Figure 3: ¹³C sensitivity (S/N) test: ¹³C spectrum of glycine (natural abundance) at 10 kHz spinning.

¹⁵N CP and sensitivity
 glycine (natural abundance), 10 kHz spinning
¹⁵N CP at 70 kHz (3.5 us) / 3 ms, plw1=82 W
 decoupling: SPINAL64 at 140 kHz (1.8 us) / 29 ms, plw12=45 W
 temperature control on: T(set) = 293 K, VT = 500 l/h

S/N: 15 (ns=64, sinocal, 20 ppm noise)




Figure 4: ¹⁵N sensitivity (S/N) test: ¹⁵N spectrum of glycine (natural abundance) at 10 kHz spinning.




Figure 5: ^1H - ^{15}N - ^{13}C double CP efficiency test: Glycine-1,2- $^{13}\text{C}_2$ - ^{15}N at 10 kHz spinning. Blue line: double CP experiment at the -1 side band CP condition for ^{13}C , red line: standard CP experiment scaled by down by a factor of 0.52




Figure 6: ^1H - ^{15}N - ^{13}C double CP efficiency test: Glycine-1,2- $^{13}\text{C}_2$ - ^{15}N at 10 kHz spinning. Blue line: double CP experiment at the +1 side band CP condition for ^{13}C , red line: standard CP experiment scaled by down by a factor of 0.44.

^1H - ^{15}N - ^{13}C double CP stability test
 glycine-1,2- ^{13}C - ^{15}N , 10 kHz spinning
 ^1H - ^{15}N CP at 70 kHz for ^{15}N , 50% to 100% ramped ^1H contact pulse, 3 ms
 ^{15}N - ^{13}C CP at 35 kHz for ^{15}N and 55 kHz for ^{13}C , square ^{15}N contact pulse, 90% to 100% ramped ^{13}C contact pulse, 20 ms
 temperature control on, T(set)=293 K, VT=500 l/h

512 spectra with 8 scans each at 5 s recycle delay




Figure 7: ^1H - ^{15}N - ^{13}C double CP stability test for glycine-1,2- ^{13}C - ^{15}N spinning at 10 kHz: 512 spectra with 8 scans each at 5 s recycle delay. Column projection (left) shows amplitude stability of CP signal, row projection (top) shows maximum noise compared to maximum signal.




Figure 8: Expanded display of the column projection in Figure 7: signal amplitude is stable over time and the noise stays constant.

¹⁵N CP stability test
 glycine (natural abundance), 10 kHz spinning
 15N CP at 70 kHz (3.5 us) / 3 ms, plw1=82 W
 decoupling: SPINAL64 at 140 kHz (1.8 us) / 29 ms, plw12=45 W
 temperature control on: T(set) = 293 K, VT = 500 l/h
 256 spectra with 4 scans each, 5.04 s repetition time




Figure 9: ¹⁵N CP stability test run on glycine (natural abundance) at 10 kHz spinning: 256 spectra with 4 scans each at 5 s recycle delay and 5.04 s repetition time. Column projection (left) shows amplitude stability of CP signal, row projection (top) shows maximum noise compared to maximum signal.




Figure 10: Expanded display of the column projection in Figure 9: signal amplitude is stable over time and the noise level stays constant.